



[www.vathin.com](http://www.vathin.com)  
[marketing@vathin.com](mailto:marketing@vathin.com)

1st floor, Building 12 , No.9, West Chuanqi Road, Jiu Hua Economic and  
Technological Development Zone, Xiangtan 411100, China

Version: RMD2502001



## Vathin Single-Use Flexible Bronchoscope Ultra-Slim Scope

Slim but mighty, ultra-slim profile to navigate  
even the narrowest passages.



### Ultra-Slim Design:

1.2mm ID, 2.6mm OD for less invasive use and easier navigation through tight airways.



### Wide Range of Applications:

Suitable for pediatric, adult, and elderly patients in diverse clinical settings, including emergency and diagnostics.




### Strong Suction Power:

High suction strength clears airways effectively, ensuring thorough examination with minimal trauma.

Ultra Slim



## Maintaining Precision By Single-use Technology



## Specification

| Model       | Working channel ID (mm) | Insertion Tube OD (mm) | Working length (mm) | Bending lock | Suction function | Bending Section | Field of view |
|-------------|-------------------------|------------------------|---------------------|--------------|------------------|-----------------|---------------|
| BC-P1600    | 3.6Fr                   | 2.6mm                  | 680                 | N            | Y                | U210°/<br>D210° | 120°±15%      |
| BC-P1601    |                         |                        | 650                 | N            |                  |                 |               |
| BC-P1602    |                         |                        | 600                 | N            |                  |                 |               |
| BC-P1603    |                         |                        | 550                 | N            |                  |                 |               |
| BC-P1604    |                         |                        | 500                 | N            |                  |                 |               |
| BC-P1600-L  |                         |                        | 680                 | Y            |                  |                 |               |
| BC-P1601-L  |                         |                        | 650                 | Y            |                  |                 |               |
| BC-P1602-L  |                         |                        | 600                 | Y            |                  |                 |               |
| BC-P1603-L  |                         |                        | 550                 | Y            |                  |                 |               |
| BC-P1604-L  |                         |                        | 500                 | Y            |                  |                 |               |
| BC-P16001   |                         |                        | 680                 | N            |                  |                 |               |
| BC-P16011   |                         |                        | 650                 | N            |                  |                 |               |
| BC-P16021   |                         |                        | 600                 | N            |                  |                 |               |
| BC-P16031   |                         |                        | 550                 | N            |                  |                 |               |
| BC-P16041   |                         |                        | 500                 | N            |                  |                 |               |
| BC-P16001-L |                         |                        | 680                 | Y            |                  |                 |               |
| BC-P16011-L |                         |                        | 650                 | Y            |                  |                 |               |
| BC-P16021-L |                         |                        | 600                 | Y            |                  |                 |               |
| BC-P16031-L |                         |                        | 550                 | Y            |                  |                 |               |
| BC-P16041-L |                         |                        | 500                 | Y            |                  |                 |               |